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HEALTHCARE®

**Kills | Tuent
C. difficile
in | en
3 minutes**

**CLEANER & DISINFECTANT / NETTOYANT et DÉSINFECTANT
BROAD SPECTRUM VIRUCIDE / VIRUCIDE À LARGE SPECTRE**

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PATHOGENS / MICROORGANISMS

- Kills 33 bacteria including:**
- Carbapenem resistant *Klebsiella pneumoniae* (CRKP)
 - Clostridium difficile*-or- *C. difficile*-or- *C. diff.* (C. diff.) spores
 - Community acquired Methicillin resistant *Staphylococcus aureus* (CA-MRSA)
 - Enterococcus faecalis*, vancomycin resistant -or- Vancomycin resistant *Enterococcus faecalis* (VRE)
 - ESBL (Extended Spectrum Beta Lactamase) producing *Klebsiella pneumoniae*
 - Escherichia coli* O157:H7
 - Escherichia coli* New Delhi Metallo-Beta Lactamase-1 (NDM-1) (NDM-1 E. coli)
 - Klebsiella pneumoniae*
 - Linezolid resistant *Staphylococcus aureus* (LRSa)
 - Listeria monocytogenes* (Listeria)
 - Legionella pneumophila*
 - Methicillin resistant *Staphylococcus aureus* (MRSA)
 - Multi-drug resistant *Enterococcus faecium* (MDR E. faecium)
 - Pseudomonas aeruginosa* (*Pseudomonas*)
 - Salmonella enterica* tested as *Salmonella choleraesuis*
 - Staphylococcus aureus*
- Kills 4 fungi including:**
- Aspergillus brasiliensis* tested as *Aspergillus niger* (mildew)
 - Trichophyton mentagrophytes* (Athlete's Foot Fungus)
- Kills 10 non-enveloped virus including:**
- Adenovirus (type 2) (Strain Adenoid 6) (common cause of respiratory infections)
 - Enterovirus type D68
 - (Human) Hepatitis A [Virus] (HAV)
 - Norovirus -or- Norwalk virus (as Feline Calicivirus)
 - Poliovirus type 1, Chat Strain
 - Rhinovirus (type 37) (a leading) causative agent of the common cold (Strain 151-1)
 - Rotavirus (causative agent of viral diarrhea)
- Kills 15 enveloped virus including:**
- Avian Influenza A (H5N1) virus
 - HIV (Human Immunodeficiency Virus) [type 1]
 - (Human) Coronavirus (additional) causative agent of the common cold
 - (Human) Hepatitis B [Virus] (HBV) [as duck HBV]
 - (Human) Hepatitis C [Virus] (HCV) [as bovine diarrhea virus]
 - (Human) Herpes virus [type 2] -or- Herpes simplex virus [type 2] (causative agent of genital herpes)
 - Influenza A [Virus] (representative of the common flu virus)
 - Measles Virus
 - MERS-Coronavirus (Middle East Respiratory Syndrome Coronavirus) [MERS-CoV]
 - Respiratory syncytial virus (RSV) (cause of respiratory infections in infants) (Strain Long)
 - SARS-Associated Coronavirus (SARS-CoV)
- Bloodborne Pathogens [virus]:**
- HIV (Human Immunodeficiency Virus) [type 1]
 - (Human) Hepatitis B [Virus] (HBV)
 - Hepatitis C [Virus] (HCV) (Bovine Viral Diarrhea virus as a surrogate virus for human Hepatitis C virus)

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- Eliminates mixing and measuring
- Strong cloth-like wipes
- Ready-to-use hospital grade disinfectant

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DIRECTIONS FOR USE: Not for cleaning or sanitizing skin. Do not use as a diaper wipe or for personal cleansing. A potable water rinse is required for food contact surfaces and for baby items and toys. Always refer to manufacturer's care instructions before using on equipment. Do not use on aluminum, brass, copper, silver, natural marble, sealed granite or chipped enamel. **DISPENSER DIRECTIONS:** Remove lid and discard seal. Insert corner of centre sheet through flap opening on lid. DO NOT push finger through opening. Replace lid and pull wipe at an angle. The next sheet pops up automatically. When finished, close lid flap to retain moisture. **To clean, disinfect and deodorize hard, nonporous surfaces:** Wipe hard, nonporous surface to be disinfected. Use enough wipes for treated surface to remain visibly wet for the contact time listed per this table. Let air dry. Gross filth must be removed prior to disinfecting. If streaking is observed, wipe with a clean, damp cloth or paper towel after appropriate contact time has expired. For surfaces and/or objects that may come into direct contact with food, a rinse with potable water is recommended. **Special Instructions for Cleaning Prior to Disinfection against *Clostridium difficile* spores:** This product is intended for use to kill or inactivate spores of *C. difficile* in settings where it is expected to be encountered. **Personal protection:** Personnel that clean items potentially soiled with spores of *C. difficile* must wear appropriate personal protective equipment (i.e., barrier protection), such as disposable gloves, gowns, masks or eye coverings. **Cleaning procedure:** Target surfaces must be cleaned prior to disinfection, as the presence of heavy soil such as fecal matter on a surface may reduce the intended disinfectant efficacy of the product. **Infectious materials disposal:** Any materials used in the cleaning process that may contain feces/waste are to be disposed of immediately in accordance with local regulations for infectious materials disposal. **Special Instructions for Cleaning Prior to Disinfection against bloodborne viral pathogens:** This product has efficacy claims against bloodborne viral pathogens (e.g., Human Immunodeficiency virus, Hepatitis B virus, Hepatitis C virus). It is intended for use against the bloodborne pathogens in settings where the viruses would be expected to be encountered. **Personal protection:** Personnel that clean items soiled with blood or body fluids must wear appropriate personal protective equipment (i.e., barrier protection), such as disposable gloves, gowns, masks or eye coverings. **Cleaning procedure:** Target surfaces must be cleaned prior to disinfection, as the presence of heavy soil on a surface may reduce the intended disinfectant efficacy of the product. **Infectious materials disposal:** Any materials used in the cleaning process that may contain blood or body fluids are to be disposed of immediately in accordance with local regulations for infectious materials disposal.

CAUTION : IRRITANT. MAY IRRITATE EYES AND SKIN. DANGEROUS FUMES FORM WHEN MIXED WITH OTHER PRODUCTS. Do not use or mix with toilet bowl cleaners, rust removers, acids or products containing ammonia. Do not get in eyes, on skin or on clothing. For sensitive skin or prolonged use, wear rubber gloves. Avoid contact with foods. **KEEP OUT OF REACH OF CHILDREN AND PETS. FIRST AID TREATMENT:** Contains sodium hypochlorite. If in eyes, rinse with water for 15 minutes. If irritation persists, call a doctor. If on skin, rinse well with water.

DISPOSAL: Dispose of wipe in trash after use. Do not flush. Offer empty container for recycling. If recycling is not available, discard container in trash. **CONTAINS:** Sodium Hypochlorite 0.55%. Contains no phosphorus. **MISE AU REBUT :** Jeter les serviettes aux ordures après l'emploi. Ne pas jeter dans les toilettes. Mettre le contenant vide au recyclage, si possible; sinon, le jeter aux ordures. **CONTIENT :** 0,55 % d'hypochlorite de sodium. Sans phosphore. **AUTHORIZED USER OF REGISTERED TRADEMARK / UTILISATEUR AUTORISÉ DE LA MARQUE DÉPOSÉE:** CLOROX HEALTHCARE. OWNED BY / QUI EST LA PROPRIÉTÉ DE THE CLOROX COMPANY, OAKLAND, CA 94612.

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Kills | Tuent
C. difficile
in | en
3 minutes

PATHOGENS / MICROORGANISMS

Kills 33 bacteria including:

- Carbapenem resistant *Klebsiella pneumoniae* (CRKP)
- Clostridium difficile* -or- *C. difficile* -or- *C. diff.* (C. diff.) spores
- Community acquired Methicillin resistant *Staphylococcus aureus* (CA-MRSA)
- Enterococcus faecalis*, vancomycin resistant -or- Vancomycin resistant *Enterococcus faecalis* (VRE)
- ESBL (Extended Spectrum Beta Lactamase) producing *Klebsiella pneumoniae*
- Escherichia coli* O157:H7
- Escherichia coli* New Delhi Metallo-Beta Lactamase-1 (NDM-1) (NDM-1 E. coli)
- Klebsiella pneumoniae*
- Linezolid resistant *Staphylococcus aureus* (LRSA)
- Listeria monocytogenes* (Listeria)
- Legionella pneumophila*
- Methicillin resistant *Staphylococcus aureus* (MRSA)
- Multi-drug resistant *Enterococcus faecium* (MDR E. faecium)
- Pseudomonas aeruginosa* (*Pseudomonas*)
- Salmonella enterica* tested as *Salmonella choleraesuis*
- Staphylococcus aureus*
- Streptococcus pneumoniae* (Strep)
- Streptococcus pyogenes* (Strep)
- Vancomycin intermediate resistant *Staphylococcus aureus* (VISA)
- Vancomycin resistant *Staphylococcus aureus* (VISA)

Mycobacterium:

- Mycobacterium bovis* (Tuberculosis or TB)
- Kills 4 fungi including:**
- Aspergillus brasiliensis* tested as *Aspergillus niger* (mildew)
- Trichophyton mentagrophytes* (Athlete's Foot Fungus)

Kills 10 non-enveloped virus including:

- Adenovirus (type 2) (Strain Adenoid 6) (common cause of respiratory infections)
- Enterovirus type D68
- (Human) Hepatitis A (Virus) (HAV)
- Norovirus -or- Norwalk virus (as Feline Calicivirus)
- Poliovirus type 1, Chat Strain
- Rhinovirus (type 37) (la leading) causative agent of the common cold (Strain 151-1)
- Rotavirus (causative agent of viral diarrhea)

Kills 15 enveloped virus including:

- Avian Influenza A (H5N1) virus
- HV (Human Immunodeficiency Virus) (type 1)
- (Human) Coronavirus (additional) causative agent of the common cold
- (Human) Hepatitis B (Virus) (HBV) (as duck HBV)
- (Human) Hepatitis C (Virus) (HCV) (as bovine diarrhea virus)
- (Human) Herpes virus (type 2) -or- Herpes simplex virus (type 2) (causative agent of genital herpes)
- Influenza A (Virus) (representative of the common flu virus)
- Measles virus
- MERS-Coronavirus (Middle East Respiratory Syndrome Coronavirus) (MERS-CoV)
- Respiratory syncytial virus (RSV) (cause of respiratory infections in infants) (Strain Long)
- SARS-Associated Coronavirus (SARS - CoV)

Bloodborne Pathogens (virus):

- HV (Human Immunodeficiency Virus) (type 1)
- (Human) Hepatitis B (Virus) (HBV)
- Hepatitis C (Virus) (HCV) (Bovine Viral Diarrhea virus as a surrogate virus for human Hepatitis C virus)

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- Eliminates mixing and measuring • Strong cloth-like wipes
- Ready-to-use hospital grade disinfectant

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DIRECTIONS FOR USE: Not for cleaning or sanitizing skin. Do not use as a diaper wipe or for personal cleansing. A potable water rinse is required for food contact surfaces and for baby items and toys. Always refer to manufacturer's care instructions before using on equipment. Do not use on aluminum, brass, copper, silver, natural marble, sealed granite, or chipped enamel. **DISPENSER DIRECTIONS:** Remove lid and discard seal. Insert corner of centre sheet through flap opening in lid. DO NOT push finger through opening. Replace lid and pull wipe at an angle. The next sheet pop up automatically. When finished, close lid flap to retain moisture. **To clean, disinfect and deodorize hard, nonporous surfaces:** Wipe hard, nonporous surface to be disinfected. Use enough wipes for treated surface to remain visibly wet for the contact time listed per this table. Let air dry. Gross fill must be removed prior to disinfecting. If streaking is observed, wipe with a clean, damp cloth or paper towel after appropriate contact time has expired. For surfaces and/or objects that may come into direct contact with food, a rinse with potable water is recommended. **Special Instructions for Cleaning Prior to Disinfection against *Clostridium difficile* spores:** This product is intended for use to kill or inactivate spores of *C. difficile* in settings where it is expected to be encountered. **Personal protection:** Personnel that clean items potentially soiled with spores of *C. difficile* must wear appropriate personal protective equipment (i.e., barrier protection), such as disposable gloves, gowns, masks or eye coverings. **Cleaning procedure:** Target surfaces must be cleaned prior to disinfection, as the presence of heavy soil such as fecal matter on a surface may reduce the intended disinfectant efficacy of the product. **Infectious materials disposal:** Any materials used in the cleaning process that may contain feces/waste are to be disposed of immediately in accordance with local regulations for infectious materials disposal. **Special Instructions for Cleaning Prior to Disinfection against bloodborne viral pathogens:** This product has efficacy claims against bloodborne viral pathogens (i.e., Human Immunodeficiency Virus, Hepatitis B virus, Hepatitis C virus). It is intended for use against the bloodborne pathogens in settings where the viruses would be expected to be encountered. **Personal protection:** Personnel that clean items soiled with blood or body fluids must wear appropriate personal protective equipment (i.e., barrier protection), such as disposable gloves, gowns, masks or eye coverings. **Cleaning procedure:** Target surfaces must be cleaned prior to disinfection, as the presence of heavy soil on a surface may reduce the intended disinfectant efficacy of the product. **Infectious materials disposal:** Any materials used in the cleaning process that may contain blood or body fluids are to be disposed of immediately in accordance with local regulations for infectious materials disposal.

CAUTION : IRRITANT. PEUT IRRITER LES YEUX ET LA PEAU. DANGEROUS FUMES FORM WHEN MIXED WITH OTHER PRODUCTS. Do not use or mix with toilet bowl cleaners, rust removers, acids, or products containing ammonia. Do not get in eyes, on skin or on clothing. For sensitive skin or prolonged use, wear rubber gloves. Avoid contact with foods. **KEEP OUT OF REACH OF CHILDREN AND PETS. FIRST AID TREATMENT:** Contains sodium hypochlorite. If in eyes, rinse with water for 15 minutes. If irritation persists, call a doctor. If on skin, rinse well with water.

DISPOSAL: Dispose of wipe in trash after use. Do not flush. Offer empty container for recycling. If recycling is not available, discard container in trash. **CONTAINS:** Sodium Hypochlorite 0.55%. Contains no phosphorus.

AUTHORIZED USER OF REGISTERED TRADEMARK / UTILISATEUR AUTORISÉ DE LA MARQUE DÉPOSÉE CLOROX HEALTHCARE, OWNED BY / OUI EST LA PROPRIÉTÉ DE THE CLOROX COMPANY, OAKLAND, CA 94612. MADE IN USA / FABRIQUÉ AUX É.-U. Imported and Distributed by / Importé et distribué par THE CLOROX COMPANY OF CANADA, LTD., 150 BISCAYNE CRES., BRAMPTON, ONTARIO L6W 4V3 Questions? 1-866-789-4973 www.cloroxhealthcare.ca

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DIRECTIONS FOR USE: Not for cleaning or sanitizing skin. Do not use as a diaper wipe or for personal cleansing. A potable water rinse is required for food contact surfaces and for baby items and toys. Always refer to manufacturer's care instructions before using an equipment. Do not use on aluminum, brass, copper, silver, natural marble, sealed granite, or chipped enamel. **DISPENSER DIRECTIONS:** Remove lid and discard seal. Insert corner of centre sheet through flap opening on lid. DO NOT push finger through opening. Replace lid and pull seal at an angle. The next sheet pops up automatically. When finished, close lid flap to retain moisture. **To clean, disinfect and deodorize hard, nonporous surfaces:** Wipe hard, nonporous surface to be disinfected. Use enough wipes for treated surface to remain visibly wet for the contact time listed per this table. Let air dry. Gross filth must be removed prior to disinfecting. If streaking is observed, wipe with a clean, damp cloth or paper towel after appropriate contact time has expired. For surfaces and/or objects that may come into direct contact with food, a rinse with potable water is recommended.

Special Instructions for Cleaning Prior to Disinfection against *Clostridium difficile* spores: This product is intended for use to kill or inactivate spores of *C. difficile* in settings where it is expected to be encountered. **Personal protection:** Personnel that clean items potentially soiled with spores of *C. difficile* must wear appropriate personal protective equipment (i.e., barrier protection), such as disposable gloves, gowns, masks or eye coverings. **Cleaning procedure:** Target surfaces must be cleaned prior to disinfection, as the presence of heavy soil such as fecal matter on a surface may reduce the intended disinfectant efficacy of the product. **Infectious materials disposal:** Any materials used in the cleaning process that may contain feces/waste are to be disposed of immediately in accordance with local regulations for infectious materials disposal. **Special Instructions for Cleaning Prior to Disinfection against bloodborne viral pathogens:** This product has efficacy claims against bloodborne viral pathogens (e.g., Human Immunodeficiency virus, Hepatitis B virus; Hepatitis C virus). It is intended for use against the bloodborne pathogens in settings where the viruses would be expected to be encountered. **Personal protection:** Personnel that clean items soiled with blood or body fluids must wear appropriate personal protective equipment (i.e., barrier protection), such as disposable gloves, gowns, masks or eye coverings. **Cleaning procedure:** Target surfaces must be cleaned prior to disinfection, as the presence of heavy soil on a surface may reduce the intended disinfectant efficacy of the product. **Infectious materials disposal:** Any materials used in the cleaning process that may contain blood or body fluids are to be disposed of immediately in accordance with local regulations for infectious materials disposal.

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CAUTION : IRRITANT. MAY IRRITATE EYES AND SKIN. DANGEROUS FUMES FORM WHEN MIXED WITH OTHER PRODUCTS. Do not use or mix with toilet bowl cleaners, rust removers, acids, or products containing ammonia. Do not get in eyes, on skin or on clothing. For sensitive skin or prolonged use, wear rubber gloves. Avoid contact with foods. **KEEP OUT OF REACH OF CHILDREN AND PETS. FIRST AID TREATMENT:** Contains sodium hypochlorite. If in eyes, rinse with water for 15 minutes. If irritation persists, call a doctor. If on skin, rinse well with water.

DISPOSAL: Dispose of wipe in trash after use. Do not flush. Offer empty container for recycling. If recycling is not available, discard container in trash. **CONTAINS:** Sodium Hypochlorite 0.55%. Contains no phosphorus. **AUTHORIZED USER OF REGISTERED TRADEMARK / UTILISATEUR AUTORISÉ DE LA MARQUE DÉPOSÉE:** CLOROX HEALTHCARE, OWNED BY / QUI EST LA PROPRIÉTÉ DE THE CLOROX COMPANY, OAKLAND, CA 94612. **Imported and Distributed by / Importé et distribué par :** THE CLOROX COMPANY OF CANADA, LTD. 150 BISCAYNE CRES., BRAMPTON, ONTARIO L6W 4V3 Questions? 1-866-789-4973 **MADE IN USA / FABRIQUÉ AUX É.-U.** www.cloroxhealthcare.ca

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SIDE SEAL AREA

SIDE SEAL AREA

Product label for Clorox Healthcare Bleach Germicidal Wipes, containing safety warnings, usage instructions, and contact information in both English and French.

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Kills bloodstream pathogens:

- Human Immunodeficiency virus type 1
- Duck Hepatitis B virus as a surrogate virus for Human Hepatitis B virus
- Bovine Viral Diarrhea virus as a surrogate virus for Human Hepatitis C virus

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DIRECTIONS FOR USE: To Clean, Deodorize and Disinfect Hard, Nonporous Surfaces in one step: Spray pull top and aim 15 cm – 20 cm from surface until surface is thoroughly wet. Allow this product to remain wet for time per the table. Then wipe with a clean, damp cloth. For heavily soiled areas, a pre-cleaning is required. Rinse with potable water for food-contact surfaces.

For Toilets: Flush toilet. Pull top and aim 15 cm – 20 cm from surface until surface is thoroughly wet. Allow this product to remain wet for the table. Brush entire bowl, including under rim, with a scrub brush or mop. Flush.

Special Instructions for Cleaning Prior to Disinfection against bloodborne viral pathogens: This product has efficacy claims against bloodborne viral pathogens (e.g., Human Immunodeficiency virus, Hepatitis B virus, Hepatitis C virus). It is intended for use against the bloodborne pathogens in settings where the viruses would be expected to be encountered, such as settings where contamination by blood or body fluids is likely. **Personal protection:** Personnel that clean items soiled with blood or body fluids must wear appropriate personal protective equipment (i.e., barrier protection), such as disposable gloves, gowns, masks or eye coverings. **Cleaning:** Target surfaces must be cleaned prior to disinfection, as the presence of heavy soil on a surface may reduce the intended disinfectant efficacy of the product. **Infectious materials disposal:** Any materials used in the cleaning process that may contain blood or body fluids are to be disposed of immediately in accordance with local regulations for infectious materials disposal.

Special Instructions for Cleaning Prior to Disinfection against Clostridium difficile spores: This product is intended for use to kill or inactivate spores of *C. difficile* in settings where it is expected to be encountered, such as settings where contamination by fecal matter is likely. **Personal protection:** Personnel that clean items potentially soiled with spores of *C. difficile* must wear appropriate personal protective equipment (i.e., barrier protection), such as disposable gloves, gowns, masks or eye coverings. **Cleaning procedure:** Target surfaces must be cleaned prior to disinfection, as the presence of heavy soil such as fecal matter on a surface may reduce the intended disinfectant efficacy of the product. **Infectious materials disposal:** Any materials used in the cleaning process that may contain feces/waste are to be disposed of immediately in accordance with local regulations for infectious materials disposal.

CAUTION: IRRITANT. MAY IRRITATE EYES AND SKIN.
Do not get in eyes, skin or clothing. For sensitive skin or prolonged use, wear rubber gloves. Avoid contamination of food. Do not mix with ammonia or acids such as vinegar, rust removers, or toilet bowl cleaners. **KEEP OUT OF REACH OF CHILDREN AND PETS. FIRST AID TREATMENT:** Contains Sodium Hypochlorite. If swallowed, call a Poison Control Centre or doctor immediately. Do not induce vomiting. If in eyes, rinse with water for 15 minutes. If on skin, rinse well with water.

STORAGE AND DISPOSAL: Store this product in a cool, dry area, away from direct sunlight and heat. **Do not refillable container.** Do not reuse or refill this container. When discarding container, place in trash collection or recycle or take to recycling center. Do not contaminate water, food, or feed by storage and disposal. **CONTAINS:** Sodium Hypochlorite 0.65%. Contains no phosphorus. **RANGEMENT ET MISE AU REBUT :** Ranger ce produit dans un endroit frais et sec, à l'abri de la chaleur et des rayons directs du soleil. Contenant non réutilisable. Ne pas réutiliser ni remplir de nouveau ce contenant vide. Mettre le contenant vide aux ordures ou au recyclage, ou à l'incinérateur à un centre de recyclage. Éviter de contaminer l'eau, ainsi que la nourriture pour les humains ou les animaux, par l'entreposage ou l'élimination. **CONTIENT :** 0,65 % d'hypochlorite de sodium. Ne contient pas de phosphore.

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AUTHORIZED USER OF REGISTERED TRADEMARK / UTILISATEUR AUTORISÉ DE LA MARQUE DÉPOSÉE CLOROX HEALTHCARE. OWNED BY / QUI EST LA PROPRIÉTÉ DE THE CLOROX COMPANY, OAKLAND, CA 94612.