Kills 44 pathogens: bacteria, virus, TB and fungi.

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**NON-BLEACH • SANS EAU DE JAVEL**

KEEP OUT  
OF REACH  
OF CHILDREN  
AND PETS.

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**NEW/OPEN VIRUSES: 30-second contact time**

- 2009-H1N1 Influenza A virus (New/2009 H1N1)
- Cytomegalovirus
- Duck Hepatitis B virus (surrogate for Human Hepatitis B virus)
- Hepatitis C virus (Bovine Viral Diarrhea virus as surrogate)
- Herpes simplex virus type 1
- Herpes simplex virus type 2
- Human Coronavirus
- Human Immunodeficiency Virus type 1 (HIV-1)
- Influenza A H5N1 virus
- Influenza A virus (2009 H1N1)
- Influenza B virus
- Respiratory Syncytial virus (RSV)

**1-minute contact time**

- Avian Influenza A virus (H5N1)
- (NON-ENVELOPED) VIRUSES: 1-minute contact time
- Rhinovirus (type 37)
- Rotavirus (Strain WA)

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**DIRECTIONS FOR USE:** This product cleans, disinfects and deodorizes hard, nonporous medical surfaces in one step with no rinsing required. When used as directed, this product contains 1.40% [14,000 ppm] of active disinfectant making it highly effective against a broad spectrum of bacteria, viruses, TB and fungi on hard, nonporous surfaces. This product is a one-step hospital-use disinfectant cleaner and deodorant designed for general cleaning and disinfecting hard, nonporous inanimate surfaces. Not for cleaning or sanitizing skin. Do not use as a diaper wipe or for personal cleansing. A potable water rinse is required for food contact surfaces. Do not use on glassware, utensils, or dishes. For surfaces and/or objects that may come into direct contact with children at the mouthing stage of development, a rinse with potable water is recommended.

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**Claims against bloodborne viral pathogens:** This product is intended for use against the bloodborne pathogens in settings where the viruses would be expected to be encountered, such as settings where contamination by blood or body fluid is likely. **Special instructions for using this product to clean and decontaminate against HIV-1, HBV and HCV on surfaces/objects soiled with blood/body fluids: Personal Protection:** When handling items soiled with blood or body fluids, use disposable impervious gloves, gowns, masks and eye coverings. **Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces and other objects before applying this product. **Contact Time:** Allow surface to remain wet for HBV, HCV and HIV-1 for 30 seconds. The contact times for other bacteria and viruses may differ. See product label for contact times. **Disposal of Infectious Materials:** Use disposable impervious gloves, gowns, masks and eye coverings. Blood and other body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal.

**CAUTION: IRRITANT. LIQUID MAY IRRITATE EYES AND SKIN.** Do not get in eyes or on clothing. For sensitive skin or prolonged use, wear rubber gloves. **KEEP OUT OF REACH OF CHILDREN AND PETS.** **FIRST AID TREATMENT:** Contains Hydrogen Peroxide. If in eyes, rinse with water for 15 minutes. Have the product container or label with you when calling a Poison Control Centre or doctor, or going for treatment.

**STORAGE AND DISPOSAL:** Store this product in a cool dry area to avoid deterioration. Dispose of wipe in trash after use. Do not flush. Recycle empty container or discard in trash.

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**U.S. Patents apply to this product. For details visit: www.thecloroxcompany.com/patents.**

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**Imported by / Importé par: THE CLOROX COMPANY OF CANADA, LTD. 150 BISCAVNE CRES., BRAMPTON, ONTARIO L6W 4V3**

**Questions? 1-866-789-4973 or / ou healthcare@clorox.com - CloroxHealthcare.ca**

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**NON-BLEACH • SANS EAU DE JAVEL**

KEEP OUT  
OF REACH  
OF CHILDREN  
AND PETS.  
DIN 02406225

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- Eliminates mixing and measuring
- For Hospital and Institutional Use
- Ready-to-Use Cleaner and Disinfectant
- Disinfection of Non-Critical Instruments

PATHOGEN	Bacteria	Virus	TB (mycobacteria)	Mold/Mildew/Fungi
CONTACT TIME	Between 30 seconds and 3 minutes	5 minutes	5 minutes	5 minutes

**DIRECTIONS FOR USE:** This product cleans, disinfects and deodorizes hard, nonporous medical surfaces in one step with no rinsing required. When used as directed, this product contains 1.40% [14,000 ppm] of active disinfectant making it highly effective against a broad spectrum of bacteria, viruses, TB and fungi on hard, nonporous surfaces. This product is a one-step hospital-use disinfectant cleaner and deodorant designed for general cleaning and disinfecting hard, nonporous inanimate surfaces. Not for cleaning or sanitizing skin. Do not use as a diaper wipe or for personal cleaning. A potable water rinse is required for food contact surfaces. Do not use on glassware, utensils, or dishes. For surfaces and/or objects that may come into direct contact with children at the mouting stage of development, a rinse with potable water is recommended.

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**Special instructions for using this product to clean and decontaminate against HIV-1, HBV and HCV on surfaces/objects soiled with blood/body fluids: Personal Protection:** When handling items soiled with blood or body fluids, use disposable impervious gloves, gowns, masks and eye coverings. **Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces and other objects before applying this product. **Contact Time:** Allow surface to remain wet for HBV, HCV and HIV-1 for 30 seconds. The contact times for other bacteria and viruses may differ. See product label for contact times. **Disposal of Infectious Materials:** Use disposable impervious gloves, gowns, masks and eye coverings. Blood and other body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal.

**CAUTION: IRRITANT. LIQUID MAY IRRITATE EYES AND SKIN.** Do not get in eyes or on clothing. For sensitive skin or prolonged use, wear rubber gloves. **KEEP OUT OF REACH OF CHILDREN AND PETS. FIRST AID TREATMENT:** Contains Hydrogen Peroxide. If in eyes, rinse with water for 15 minutes. Have the product container or label with you when calling a Poison Control Centre or doctor, or going for treatment.

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# Hydrogen Peroxide Cleaner Disinfectant Wipes

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KEEP OUT OF REACH OF CHILDREN AND PETS.

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**DIRECTIONS FOR USE:** This product cleans, disinfects and deodorizes hard, nonporous medical surfaces in one step with no rinsing required. When used as directed, this product contains 1.40% [14,000 ppm] of active disinfectant making it highly effective against a broad spectrum of bacteria, viruses, TB and fungi on hard, nonporous surfaces. This product is a one-step hospital-use disinfectant cleaner and deodorant designed for general cleaning and disinfecting hard, nonporous inanimate surfaces. Not for cleaning or sanitizing skin. Do not use as a diaper wipe or for personal cleansing. A potable water rinse is required for food contact surfaces. Do not use on glassware, utensils, or dishes. For surfaces and/or objects that may come into direct contact with children at the mounting stage of development, a rinse with potable water is recommended.

**To Clean and Disinfect Hard, Nonporous Surfaces:** Wipe surface until completely wet. Allow this product to remain wet for contact time on label. Let air dry. For heavily soiled areas, a precleaning is required.

**Claims against bloodborne viral pathogens:** This product is intended for use against the bloodborne pathogens in settings where the viruses would be expected to be encountered, such as settings where contamination by blood or body fluid is likely.

**Special instructions for using this product to clean and decontaminate against HIV-1, HBV and HCV on surfaces/objects soiled with blood/body fluids:**

**Personal Protection:** When handling items soiled with blood or body fluids, use disposable impervious gloves, gowns, masks and eye coverings. **Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces and other objects before applying this product. **Contact Time:** Allow surface to remain wet for HBV, HCV and HIV-1 for 30 seconds. The contact times for other bacteria and viruses may differ. See product label for contact times. **Disposal of Infectious Materials:** Use disposable impervious gloves, gowns, masks and eye coverings. Blood and other body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal.

**CAUTION: IRRITANT. LIQUID MAY IRRITATE EYES AND SKIN.** Do not get in eyes or on clothing. For sensitive skin or prolonged use, wear rubber gloves. **KEEP OUT OF REACH OF CHILDREN AND PETS.**

**FIRST AID TREATMENT:** Contains Hydrogen Peroxide. If in eyes, rinse with water for 15 minutes. Have the product container or label with you when calling a Poison Control Centre or doctor, or going for treatment.

**STORAGE AND DISPOSAL:** Store this product in a cool dry area to avoid deterioration. Dispose of wipe in trash after use. Do not flush. Nonrefillable container. Do not reuse or refill empty container. Recycle empty container or discard in trash.

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HEALTHCARE®

# Hydrogen Peroxide Cleaner Disinfectant

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**KILLS 33 PATHOGENS IN 30 SECONDS**  
**TUE 33 AGENTS PATHOGENES EN 30 SECONDES**

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**GUARANTEE / GARANTIE**

Hydrogen peroxide / Peroxyde d'hydrogène ..... 1.4%

**WARNING: CONTAINS THE ALLERGEN SULFITES.**

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Contains salicylic acid at 0.1% as a preservative.

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**DIN 02403528**

REGISTRATION NO. 32403 PEST CONTROL PRODUCTS ACT

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HINGED AREA, NO PRINT ALLOWED

### DIRECTIONS FOR USE:

This product is a one-step hospital-use disinfectant, cleaner and deodorant for hard, nonporous inanimate surfaces and is also a soft surface sanitizer. When used as directed, this product contains 1.4% [14,000 ppm] of Hydrogen Peroxide making it highly effective against a broad spectrum of bacteria, viruses, TB and fungi on hard and soft nonporous surfaces.

#### To clean and disinfect hard, nonporous surfaces:

Spray on surface to remain wet for contact times listed on the label. Then wipe. For heavily soiled areas, a precleaning is required. Rinse with potable water for food-contact surfaces. For all others, no rinsing is required.

#### To clean and deodorize toilets:

To clean and deodorize toilet bowl, squirt liberally on toilet sides and upper toilet bowl rim. Brush all surfaces and flush.

**KILLS HIV, HBV AND HCV ON PRE-CLEANED, NONPOROUS SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV), and Hepatitis C Virus (HCV).

#### Special instructions for cleaning and decontamination against HIV-1 (AIDS Virus), HBV and HCV on surfaces/objects soiled with blood/body fluids:

**Personal Protection:** When handling items soiled with blood or body fluids, use disposable gloves, gowns, masks, and eye coverings.

#### Cleaning Procedure:

Blood and other body fluids must be thoroughly cleaned from surfaces and objects before applying this product.

#### Contact Time:

Let stand for 30 seconds. For all other organisms, see directions for contact time. See product label for contact times.

#### Disposal of infectious materials:

Use disposable impervious gloves, gowns, masks and eye coverings. Blood and other body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal.

#### To sanitize soft surfaces:

Hold bottle upright 15-20 cm from surface. Spray until fabric is wet. DO NOT SATURATE. Let stand for 30 seconds. Allow to air dry. For difficult odours or heavy fabrics, repeat application. Heavy soil must be removed prior to sanitizing.

**STORAGE:** Store this product in a cool, dry area, away from direct sunlight and heat.

### 44 PATHOGENS / MICROORGANISMS

#### BACTERIA:

##### 30-second contact time

- *Acinetobacter baumannii* (*A. baumannii*)
- *Campylobacter jejuni* (*C. jejuni*)
- Carbapenem resistant *Klebsiella pneumoniae*
- *Enterobacter aerogenes*
- *Escherichia coli* (*E. coli*)
- *Escherichia coli* O157:H7 (*E. coli*)
- *Escherichia coli* New Delhi Metallo-Beta Lactamase-1 (NDM-1)
- *Escherichia coli* with extended beta-lactamase resistance
- *Klebsiella oxytoca*

- *Klebsiella pneumoniae* New Delhi Metallo-Beta Lactamase-1 (NDM-1)
- *Klebsiella pneumoniae* with extended beta-lactamase resistance
- Methicillin resistant *Staphylococcus aureus* (Hospital acquired MRSA)
- Methicillin resistant *Staphylococcus aureus* (MRSA)

- Multi-drug resistant *Enterococcus faecium*

- Multi-drug resistant *Klebsiella pneumoniae*
- Penicillin resistant *Streptococcus pneumoniae* (PRSP)
- *Pseudomonas aeruginosa* (*P. aeruginosa*)
- *Salmonella enterica* (*S. enterica*)
- *Stenotrophomonas maltophilia* (*S. maltophilia*)
- *Streptococcus pyogenes* (*S. pyogenes*)
- Vancomycin resistant *Enterococcus faecalis* (VRE)

- 45-second contact time
- Community acquired Methicillin resistant *Staphylococcus aureus*
- Linezolid resistant *Staphylococcus aureus*
- *Staphylococcus epidermidis*
- 1-minute contact time
- *Staphylococcus aureus*

- MYCOBACTERIA:
- 4-minute contact time
- *Mycobacterium bovis* (BCG or TB)

#### VIRUSES:

##### 30-second contact time

- H1N1 Influenza A virus (2009 - Novel H1N1)
- Avian Influenza A virus (H5N1)
- Cytomegalovirus
- Duck Hepatitis B virus (surrogate for Human Hepatitis B virus)
- Hepatitis C virus (Bovine Viral Diarrhea virus as surrogate)
- *Herpes simplex virus type 1*
- *Herpes simplex virus type 2*
- Human Coronavirus
- Human Immunodeficiency Virus type 1 (HIV-1)
- Influenza A (H3N2) virus
- Influenza B virus
- Respiratory Syncytial virus (RSV)

- 1-minute contact time
- Feline Calicivirus
- Norovirus
- Rhinovirus (Type 37)
- Rotavirus

- MOLD, FUNGI:
- 3-minute contact time
- *Candida albicans*
- *Trichophyton mentagrophytes*

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- **4-minute contact time**
- *Mycobacterium bovis* (BCG or TB)

**VRUSES:**

- **30-second contact time**
- H1N1 Influenza A virus (2009 - Novel H1N1)
- Avian Influenza A virus (H5N1)
- Dytomégalo-virus
- Duck Hepatitis B virus (surrogate for Human Hepatitis B virus)
- Hepatitis C virus (Bovine Viral Diarrhea virus as surrogate)
- Herpes simplex virus type 1
- Herpes simplex virus type 2
- Human Coronavirus
- Human Immunodeficiency Virus type 1 (HIV-1)
- Influenza A (H3N2) virus
- Influenza B virus
- Respiratory Syncytial virus (RSV)
- **1-minute contact time**
- Feline Calicivirus
- Norovirus
- Rhinovirus (Type 37)
- Rotavirus

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**GUARANTEE / GARANTIE**

Hydrogen Peroxide / Peroxyde d'hydrogène ..... 1.4%

WARNING: CONTAINS THE ALLERGEN SULFITES.

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**DIRECTIONS FOR USE:** This product is a one-step hospital-use disinfectant, cleaner and deodorant for hard, nonporous inanimate surfaces and is also a soft surface sanitizer. When used as directed, this product contains 1.4% (14,000 ppm) of Hydrogen Peroxide making it highly effective against a broad spectrum of bacteria, viruses, TB and fungi on soft, hard, nonporous surfaces. To clean and disinfect hard, nonporous surfaces: Spray on surface to remain wet for contact times listed on the label. Then wipe. For heavily soiled areas, a pre-cleaning is required. Rinse with potable water for food-contact surfaces. For all others, no rinsing is required. To clean and deodorize toilets: To clean and deodorize toilet bowls, squirt liberally on toilet sides and upper toilet bowl rim. Brush all surfaces and flush. **KILLS HIV, HBV AND HCV ON PRE-CLEANED, NONPOROUS SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS** in healthcare settings or other settings in which there is an expected likelihood of soiling of surfaces/objects with blood/body fluids, and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential, for transmission of Human Immunodeficiency Virus Type 1 (HIV-1) (associated with AIDS), Hepatitis B Virus (HBV), and Hepatitis C Virus (HCV). **Special Instructions for cleaning and decontamination against HIV-1 (AIDS Virus), HBV and HCV on surfaces/objects soiled with blood/body fluids: Personal Protection:** When handling items soiled with blood or body fluids, use disposable gloves, gowns, masks, and eye coverings. **Cleaning Procedure:** Blood and other body fluids must be thoroughly cleaned from surfaces and objects before applying this product. **Contact Time:** Let stand for 30 seconds. For all other organisms, see directions for contact time. See product label for contact times. **Disposal of infectious materials:** Use disposable impervious gloves, gowns, masks and eye coverings. Blood and other body fluids must be autoclaved and disposed of according to local regulations for infectious waste disposal. To sanitize soft surfaces: Hold bottle upright 15-20 cm from surface. Spray until fabric is wet. **DO NOT SATURATE.** Let stand for 30 seconds. Allow to air dry. For difficult odours or heavy fabrics, repeat application. Heavy soil must be removed prior to sanitizing. **STORAGE:** Store this product in a cool, dry area, away from direct sunlight and heat.

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U.S. patents apply to this product. For details visit: www.thecloroxcompany.com/patents.

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THE CLOROX COMPANY OF CANADA, LTD., 150 BISCAYNE CRES., BRAMPTON, ONTARIO L6W 4V3

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CLOROX HEALTHCARE, OWNED BY / QUI EST LA PROPRIÉTÉ DE THE CLOROX COMPANY, OAKLAND, CA 94612.

Questions? 1-866-789-4973 or / ou healthcare@clorox.com - CloroxHealthcare.ca

based on  
**aseptix**  
 technology

135817.030

**POSITION ONLY**  
**85%**  
**0 55500 01664 8**  
 00000 00000

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**AUTHORIZED USER OF REGISTERED TRADEMARK / UTILISATEUR AUTORISÉ DE LA MARQUE DÉPOSÉE CLOROX HEALTHCARE, OWNED BY / QUI EST LA PROPRIÉTÉ DE THE CLOROX COMPANY, OAKLAND, CA 94612.**
AUTHORIZED USER OF REGISTERED TRADEMARK / UTILISATEUR AUTORISÉ DE LA MARQUE DÉPOSÉE CLOROX HEALTHCARE, OWNED BY / QUI EST LA PROPRIÉTÉ DE LA CLOROX COMPANY, OAKLAND, CA 94612.
AUTHORIZED USER OF REGISTERED TRADEMARK / UTILISATEUR AUTORISÉ DE LA MARQUE
AUTHORIZED USER OF REGISTERED TRADEMARK / UTILISATEUR AUTORISÉ DE LA MARQUE DÉPOSÉE