











Clorox Healthcare Equipment Compability Guide



Vital Sign Monitors/ Consoles

		Bleach					Hydrogen Peroxide	Quat Alcohol
								
		Clorox Healthcare® Bleach Germicidal Wipes	Clorox Healthcare® Bleach Germicidal Cleaners	Clorox Healthcare® Fuzion Cleaner Disinfectant	Dispatch® Hospital Cleaner Disinfectant Towels with Bleach	Clorox® Germicidal Bleach	Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes	Clorox Healthcare® Multi-Surface Quat Alcohol Cleaner Disinfectant Wipes
Manufacturer	Product							
Barco	Eonis® Diagnostic Displays							
General Electric Company	Transport Pro							
Midmark	IQvitals							
	IQvitals® Blood Pressure Cuffs							
Phillips	IntelliVue Patient Mx400/450/500/600/700/800							
	IntelliVue Telemetry Transmitters							
	IntelliVue Monitors							
	Avalon Fetal Monitors							
Welch Allyn	Connex Vital Signs							
	Accessory Cable Management Mobile Stand 4800/4900 Series							
	Propaq Encore® Blood Pressure Cuffs							
	Harvey™ Elite® Stethoscopes							
	LumiView™ Binocular Microscope							
	Digital MacroView™ Oscopes and Halogen HPX® Illuminators							
MindRay	Passport Series/Accutorr/Telemetry							

Legend	
	Product is mentioned by name in Cleaning & Care guide and is approved for use
	Active ingredient/level is listed in Cleaning & Care guide as approved for use, but product is not listed by name
	* Product mentioned in company's Cleaning & Care guide has the same EPA sub-registration